

---

# UNIT 9 APPLICATIONS OF EMF MEASUREMENTS

---

## Structure

- 9.1 Introduction
  - Objectives
- 9.2 Experiment 8 : Determination of pH by EMF Method
  - Principle
  - Requirements
  - Procedure
  - Observations
  - Calculations
  - Result
- 9.3 Experiment 9 : Determination of Solubility Product of AgCl from EMF Measurements
  - Principle
  - Requirements
  - Procedure
  - Observations
  - Calculations
  - Result
- 9.4 Summary

---

## 9.1 INTRODUCTION

---

In Unit 3 you were introduced to elementary concepts of potentiometry such as types of electrodes, electrode potential, galvanic cell, cell potential, emf, Nernst equation, etc. We also described a demonstration experiment for construction of a Daniell cell and measurement of its emf. In this unit we shall describe two experiments which are based on the applications of emf measurements. In these experiments, you will measure the emf of appropriate galvanic cells using low cost potentiometer and then determine the pH of solutions and the solubility product of a sparingly soluble salt such as silver chloride.

After studying this unit and having performed the experiments described in this unit, you should be able to:

- explain the terms electrodes, electrode potential, galvanic cell, cell potential, emf, Nernst equation, pH and solubility product,
- state the relationship between emf of a cell (or half-cell potential) and the concentration of species involved in cell reaction (or half-cell reaction),
- describe the relationship between emf of a cell and the pH of the solution in the indicator half-cell,
- construct copper, silver and quinhydrone electrodes,
- measure emf of a cell using low cost potentiometer,
- determine pH of an acid solution and solubility product of a sparingly soluble salt by measuring emf with a low cost potentiometer.

---


## 9.2 EXPERIMENT 8 : DETERMINATION OF pH BY EMF METHOD

---

In Chemistry Lab - I, CHE-03(L), course, you performed an experiment based on low cost pH meter. You performed a pH metric titration and determined the amount of acetic acid in a commercial sample of vinegar. The low cost pH meter can also be used for measuring pH of solutions directly. But in this experiment, you will use the low cost pH meter as a potentiometer. You will measure the emf of a galvanic cell and then, calculate the pH of hydrochloric acid and acetic acid solution by making use of Nernst equation.

## 9.2.1 Principle

In Unit 3 of this course, you studied that the quinhydrone electrode is reversible with respect to the  $H^+$  ions. You have also studied that the half-cell reaction of this electrode is



The expression for electrode potential of a quinhydrone electrode ( $E_{Q/H_2Q}$ ) in a solution of  $H^+$  ions can be written, by using Nernst equation (Eq.3.18), as follows :

$$E_{Q/H_2Q} = E_{Q/H_2Q}^0 + \frac{2.303 RT}{2F} \log [H^+]^2 \quad \dots (3.18)$$


$$= E_{Q/H_2Q}^0 + \frac{2.303 RT \times 2}{2F} \log [H^+]$$

$$= E_{Q/H_2Q}^0 + \frac{2.303 RT}{F} \text{pH} \quad \dots (3.30)$$

The standard reduction potential,  $E_{Q/H_2Q}^0$ , for quinhydrone electrode is 0.699 V. Therefore, from the above we get,

$$E_{Q/H_2Q} = 0.699 - \frac{2.303 RT}{F} \text{pH} \quad \dots (9.1)$$

Thus, you can see that the electrode potential of a quinhydrone electrode is related to the pH of the solution with which it is in contact. Now if you set up a galvanic cell by connecting the quinhydrone indicator electrode to a saturated silver-silver chloride reference electrode, the quinhydrone electrode will act as a reduction electrode and the silver-silver chloride electrode will act as an oxidation electrode. The cell can be represented as follows:



The half-cell reactions are :



Using Eq. 3.14 of Unit 3 of this course, the emf of the cell can be written as:

$$E_{\text{cell}} = E_{Q/H_2Q} - E_{AgCl/Cl^-}$$

But  $E_{AgCl/Cl^-}$  (as per literature value) = 0.222 V

Using this in Eq. 9.1, we get

$$E_{\text{cell}} = 0.699 - \frac{2.303 RT}{F} \text{pH} - 0.222$$

$$\text{or } E_{\text{cell}} = 0.477 - \frac{2.303 RT}{F} \text{pH}$$

$$\text{or } \text{pH} = \frac{(E_{\text{cell}} - 0.477)F}{2.303 RT} \quad \dots (9.2)$$

Thus, by measuring the emf of the cell, the pH of the solution can be calculated.

## 9.2.2 Requirements

Apparatus		Chemicals
1.	Low cost potentiometer	1 No.
2.	Silver-silver chloride electrode	1 No.
3.	Carbon electrode	1 No.
		1. Quinhydrone
		2. Potassium chloride

4. Salt Bridge 1 No.  
 5. Beakers 100 cm<sup>3</sup> 2 Nos.

### Stock Solutions Provided

1. Hydrochloric acid ( $1.00 \times 10^{-2}$  M)  
 2. Acetic acid - ( $1.00 \times 10^{-2}$  M)

### 9.2.3 Procedure

Calibrate the low cost potentiometer by following the method given in Sub-section 3.4.6 of Block 1 of this course.

Set up a silver-silver chloride half-cell by taking 20-25 cm<sup>3</sup> of saturated solution of KCl in one 100 cm<sup>3</sup> beaker and by placing a silver-silver chloride electrode in it. Set up a quinhydrone half-cell by taking 20-25 cm<sup>3</sup> of  $1.00 \times 10^{-2}$  M hydrochloric acid in another 100 cm<sup>3</sup> beaker. Add a pinch of quinhydrone and insert a carbon electrode in this solution. Connect the two half-cells through a KCl salt bridge. If a salt bridge is not available, take a strip of filter paper, dip it in saturated KCl solution and use it in place of the salt bridge.

Connect the carbon electrode to the indicator terminal (through the red crocodile clip) and the silver-silver chloride electrode to the reference terminal (through the black crocodile clip) of the instrument.

Keep the range selector at 2 V and then measure the emf of the cell. Remember each division on the meter scale is equal to 0.04 V in the 2 V range. Take the measurements at least two more times for the same set-up, and record the values in Table 9.1.

Now set up the quinhydrone half-cell again, this time using acetic acid ( $1.00 \times 10^{-2}$  M), and measure the emf value in the same manner. Repeat the measurement at least two more times and record the values in Table 9.1.

### 9.2.4 Observations

Room temperature,  $T = \quad ^\circ\text{C} = \quad \text{K}$

Table 9.1 : Emf Values of Cells

Cell	Emf/V
$\text{Ag}(s)   \text{AgCl}(s), \text{KCl}(\text{satd})     \text{Q}, \text{H}_2\text{Q}(s), \text{HCl}(1.00 \times 10^{-2} \text{ M})   \text{C}(s)$	----- ----- ----- -----
Average value of emf of the cell	-----
$\text{Ag}(s)   \text{AgCl}(s), \text{KCl}(\text{satd})     \text{Q}, \text{H}_2\text{Q}(s), \text{CH}_3\text{COOH}(1.00 \times 10^{-2} \text{ M})   \text{C}(s)$	----- ----- ----- -----
Average value of emf of the cell	-----

### 9.2.5 Calculations

pH value of  $1.00 \times 10^{-2}$  M HCl

Average value of emf of the cell  
 containing  $1.00 \times 10^{-2}$  M HCl = ..... V

Temperature ( $T$ ) = ..... K

Gas constant ( $R$ ) =  $8.314 \text{ J K}^{-1} \text{ mol}^{-1}$

Faraday constant ( $F$ ) =  $96500 \text{ C mol}^{-1}$

$$\text{Using Eq. 9.2, pH of } 1.00 \times 10^{-2} \text{ M HCl} = \frac{(E_{\text{cell}} - 0.477) 96500}{2.303 \times 8.314 T}$$

(You can compare this with the theoretical pH value of 2)

### pH Value of $1.00 \times 10^{-2} \text{ M CH}_3\text{COOH}$

Average value of emf of the cell containing  $1.00 \times 10^{-2} \text{ M CH}_3\text{COOH} = \dots \text{ V}$

$$\text{Again, using Eq. 9.2, pH of } 1.00 \times 10^{-2} \text{ M CH}_3\text{COOH} = \frac{(E_{\text{cell}} - 0.477) 96500}{2.303 \times 8.314 T}$$

### 9.2.6 Result

The pH values of  $1.00 \times 10^{-2} \text{ M HCl}$  and  $1.00 \times 10^{-2} \text{ M CH}_3\text{COOH}$  are found to be --- and ---, respectively.

Give reason for the relative magnitudes of pH of  $1.00 \times 10^{-2} \text{ M HCl}$  and  $1.00 \times 10^{-2} \text{ M CH}_3\text{COOH}$  in the space given below:

.....

.....

.....

.....


.....

## 9.3 EXPERIMENT 9 : DETERMINATION OF SOLUBILITY PRODUCT OF AgCl FROM EMF MEASUREMENTS

In Experiment 8 you have made use of emf of a cell in determining the pH of a solution. In Unit 17 of the CHE-04 course, you had studied that the emf measurement can also be used in determining the solubility product of a sparingly soluble salt. In this experiment, you will determine the solubility product of a sparingly soluble salt, silver chloride, by measuring the emf of an appropriate galvanic cell.

### 9.3.1 Principle

You know that the silver chloride with a solubility of  $1.37 \times 10^{-5} \text{ mol dm}^{-3}$  in water is a sparingly soluble salt. As a result, when silver chloride is placed in water, only a very small amount of it dissolves in water to give a saturated solution. You also know that all the sparingly soluble salts are strong electrolytes; therefore, AgCl that dissolves will be completely dissociated into  $\text{Ag}^+$  and  $\text{Cl}^-$  ions. Thus, there will be an equilibrium between undissolved AgCl and its ions in solution:




The equilibrium constant for this reaction, which is generally known as solubility product constant,  $K_{\text{sp}}$ , is given by

$$K_{\text{sp}} = [\text{Ag}^+(aq)][\text{Cl}^-(aq)] \quad \dots (9.3)$$

where  $[\text{Ag}^+(aq)]$  and  $[\text{Cl}^-(aq)]$  represent the molar concentrations of respective ions. The potentiometric determination of  $K_{\text{sp}}$  is based on the following simple idea.

An aqueous solution of AgCl contains equal amounts of  $\text{Ag}^+$  and  $\text{Cl}^-$  ions. However, if to this solution, a known excess of KCl is added, then  $[\text{Cl}^-]$  is essentially determined by the

concentration of KCl. The  $[Ag^+]$  can be determined by setting up and then measuring the emf of the following galvanic cell:



The emf of the cell is given by


$$\begin{aligned} E_{cell} &= E_{Ag^+/Ag} - E_{Cu^{2+}/Cu} \\ &= \left( E^\circ_{Ag^+/Ag} + \frac{2.303 RT}{F} \log [Ag^+] \right) - \left( E^\circ_{Cu^{2+}/Cu} + \frac{2.303 RT}{F} \log [Cu^{2+}] \right) \\ &= E^\circ_{Ag^+/Ag} - E^\circ_{Cu^{2+}/Cu} - \frac{2.303 RT}{2F} \log(1.000 \times 10^{-1}) + \frac{2.303 RT}{F} \log [Ag^+] \\ &= \left( E^\circ_{Ag^+/Ag} - E^\circ_{Cu^{2+}/Cu} + \frac{2.303 RT}{2F} \right) + \frac{2.303 RT}{F} \log [Ag^+] \dots (9.5) \end{aligned}$$

$$\begin{aligned} & - \frac{2.303 RT}{2F} \log(1.000 \times 10^{-1}) \\ &= - \frac{2.303 RT}{2F} (-1) \\ &= \frac{2.303 RT}{2F} \end{aligned}$$

At a given temperature, the value of the term within brackets is constant. Thus, by measuring the emf of the cell, the value of  $[Ag^+]$  can be calculated using Eq. 9.5. A more accurate method of obtaining  $[Ag^+]$  in solutions of unknown concentration from  $E_{cell}$  values is to plot a calibration curve between  $[Ag^+]$  values of silver nitrate solutions of known concentration and the corresponding  $E_{cell}$  values. This is explained below.

### Plotting a calibration curve

For the purpose of drawing calibration curve, you shall use the following cell:



where  $Ag^+(c)$  is the concentration of standard silver nitrate solution. You shall use  $AgNO_3$  solutions of  $5.000 \times 10^{-2} M$ ,  $1.000 \times 10^{-2} M$ ,  $5.000 \times 10^{-3} M$ ,  $1.000 \times 10^{-3} M$ ,  $5.000 \times 10^{-4} M$  and  $1.000 \times 10^{-4} M$  concentration and measure  $E_{cell}$  every time and note down the values in Table 9.2.

The actual plot as per Eq. 9.5 should be between  $E_{cell}$  and  $\log [Ag^+]$ . Since  $[Ag^+]$  in various  $AgNO_3$  solutions is much less than 1,  $\log [Ag^+]$  will be a negative quantity in each case. But it will be easier to plot a graph using positive quantities. Hence, Eq. 9.5 is transformed into Eq. 9.7 by multiplying the last term in right hand side of Eq. 9.5 by  $(-1)$  twice.

$$E_{cell} = \left( E^\circ_{Ag^+/Ag} - E^\circ_{Cu^{2+}/Cu} + \frac{2.303 RT}{2F} \right) - \frac{2.303 RT}{F} (-\log [Ag^+]) \dots (9.7)$$

You should use this equation for the purpose of obtaining the calibration curve. Using  $E_{cell}$  and  $-\log [Ag^+]$  values from Table 9.2, plot the calibration curve. Shape of the curve should be as shown in Fig. 9.1.




Fig. 9.1 : Calibration curve and  $[Ag^+]$  by interpolation

You can note that Eq. 9.5 is applicable to both the cells set up as per Eq. 9.4 and 9.6 since  
i) the reference electrode is the same in each case, and  
ii) the indicator electrode in each case is reversible with respect to  $Ag^+$  ions.

To understand this clearly, you are advised to go through Section 17.8 of Unit 17 of the CHE-04 course.

### Obtaining unknown $[Ag^+]$ from the calibration curve

The emf of the cell set up as per Eq. 9.4 is measured to obtain  $[Ag^+]$  in a saturated solution of AgCl by interpolation method. You have studied the details of interpolation method in Section 2.6 of Unit 2 of this course. Following this method, find out the  $x$ -coordinate value ( $b$ ) from the curve corresponding to the  $E_{cell}$  value ( $a$ ) obtained for the cell mentioned in Eq. 9.4. Once the  $x$ -coordinate value ( $b$ ) is known,  $[Ag^+]$  in a saturated solution of AgCl can be calculated using the following relationship:

$$-\log [Ag^+] = b$$

$$\text{Hence, } [Ag^+] = 10^{-b} \text{ M} \quad \dots (9.8)$$

(as shown in Example 8 discussed in Sub-section 2.3.2 of Unit 2 of Block 1 of this course).

### Calculation of $K_{sp}$

$[Cl^-] = 3.333 \times 10^{-2} \text{ M}$  as shown in the calculation part of this experiment.

Using  $[Ag^+]$  value obtained above by interpolation method and  $[Cl^-]$  as  $3.333 \times 10^{-2} \text{ M}$ , we can calculate  $K_{sp}$  using Eq. 9.3.

$$K_{sp} = ([Ag^+] \text{ as obtained by interpolation method}) \times ([Cl^-] \text{ in the indicator half-cell})$$

$$\text{or } K_{sp} = (10^{-b} \times 3.333 \times 10^{-2}) \quad \dots (9.9)$$

### 9.3.2 Requirements

#### Apparatus


1.	Low cost potentiometer	1 No.
2.	Copper electrode	1 No.
3.	Silver electrode	1 No.
4.	Salt bridge/ $KNO_3$	1 No.
5.	Beaker $100 \text{ cm}^3$	4 No.
6.	Graduated pipette $20 \text{ cm}^3$	1 No.
7.	Graduated pipette $10 \text{ cm}^3$	1 No.

#### Stock Solutions Provided

- $1.000 \times 10^{-4} \text{ M}$ ,  $5.000 \times 10^{-4} \text{ M}$ ,  $1.000 \times 10^{-3} \text{ M}$ ,  $5.000 \times 10^{-3} \text{ M}$ ,  $1.000 \times 10^{-2} \text{ M}$ ,  $5.000 \times 10^{-2} \text{ M}$  and  $1.000 \times 10^{-1} \text{ M}$   $AgNO_3$  solutions
- $1.000 \times 10^{-1} \text{ M}$   $CuSO_4$  solution
- $1.000 \times 10^{-1} \text{ M}$   $KCl$  solution

### 9.3.3 Procedure

- Calibrate the low cost potentiometer as you did in Experiment 8 by following the method given in Sub-section 3.4.6 of Block 1 of this course.
- Set up the following galvanic cell (as per Eq. 9.6) using  $1.000 \times 10^{-1} \text{ M}$   $CuSO_4$  and  $1.000 \times 10^{-4} \text{ M}$   $AgNO_3$  solutions, copper and silver electrodes, and  $KNO_3$  salt bridge:



- Connect the copper electrode to the reference terminal (through the black crocodile clip) and the silver electrode to the indicator terminal (through the red crocodile clip) of the instrument.
- Keep the range selector at 2 V and then measure the emf of the cell. Remember each division on the meter scale is equal to 0.04 V in the 2 V range. Note down the measurement in Table 9.2 and remove the salt bridge as soon as possible.

5. Now take  $5.000 \times 10^{-4}$  M silver nitrate solution in the indicator half-cell. Connect the two half-cells through the salt bridge and note down the emf. Similarly measure the emf for the rest of the silver nitrate solutions ( $1.000 \times 10^{-3}$  M,  $5.000 \times 10^{-3}$  M,  $1.000 \times 10^{-2}$  M and  $5.000 \times 10^{-2}$  M) and record the values in Table 9.2.
6. Prepare a system by mixing  $10 \text{ cm}^3$  of  $1.000 \times 10^{-1}$  M  $\text{AgNO}_3$  solution with  $20 \text{ cm}^3$  of  $1.000 \times 10^{-1}$  M KCl solution. Shake the mixture well. Take this mixture in the indicator half-cell and set up a galvanic cell as shown in Eq. 9.4.
7. Measure the emf of the cell and record the value below Table 9.2.
8. Repeat the steps 6 and 7 once more and record the value of emf.

### 9.3.4 Observations

Room temperature = °C

Table 9.2 : Emf Values of Cells constructed as per Eq. 9.8

S.No.	$[\text{Ag}^+]/\text{M}$	$E_{\text{cell}}/\text{V}$	$\log [\text{Ag}^+]/\text{M}$	$-\log [\text{Ag}^+]/\text{M}$
1.	$1.000 \times 10^{-4}$			
2.	$5.000 \times 10^{-4}$			
3.	$1.000 \times 10^{-3}$			
4.	$5.000 \times 10^{-3}$			
5.	$1.000 \times 10^{-2}$			
6.	$5.000 \times 10^{-2}$			

M stands for molarity unit

Emf of the cell constructed as per Eq. 9.4 is measured twice and the values are recorded below.

- (i) ..... V  
 (ii) ..... V

The average of the two measurements = ..... V

We shall use this as  $a$ , for the purpose of interpolation.

### 9.3.5 Calculations

#### (i) Plotting of calibration curve

Plot a calibration curve between  $E_{\text{cell}}$  and  $-\log [\text{Ag}^+]$  as explained in Sub-section 9.3.1. For this purpose, use the values from Table 9.2. You will get a calibration curve as shown in Fig. 9.1.

#### (ii) Obtaining $[\text{Ag}^+]$ from the calibration curve

As explained in Sub-section 9.3.1, obtain  $x$ -coordinate ( $b$ ) corresponding to the average  $E_{\text{cell}}$  value ( $a$ ) of the cell constructed as per Eq. 9.4.

$b = \dots\dots\dots$

Then using Eq. 9.8,  $[\text{Ag}^+] = 10^{-b}\text{M}$   
 $= \dots\dots\dots \text{M}$

(iii) Calculation of  $[Cl^-]$

$$\begin{aligned} \text{Volume of } 1.000 \times 10^{-1} \text{ M KCl taken} \\ \text{in the indicator half-cell} &= 20 \text{ cm}^3 \end{aligned}$$

$$\text{Amount of KCl present} = \frac{1.000 \times 10^{-1} \times 20}{1000} = 2.000 \times 10^{-3} \text{ mol}$$

$$\text{Volume of } 1.000 \times 10^{-1} \text{ M AgNO}_3 \text{ taken} = 10 \text{ cm}^3$$

$$\text{Amount of AgNO}_3 \text{ present} = 1.000 \times 10^{-1} \times \frac{10}{1000} = 1.000 \times 10^{-3} \text{ mol}$$

Now,  $1.000 \times 10^{-3}$  mol of KCl will react with  $1.000 \times 10^{-3}$  mol of AgNO<sub>3</sub> to give  $1.000 \times 10^{-3}$  mol of AgCl and  $1.000 \times 10^{-3}$  mol of KNO<sub>3</sub>.

Hence, amount of KCl left in the solution

$$= 2.000 \times 10^{-3} \text{ mol} - 1.000 \times 10^{-3} \text{ mol}$$

$$= (2.000 - 1.000) \times 10^{-3} \text{ mol}$$

$$= 1.000 \times 10^{-3} \text{ mol}$$

$$\text{Volume of the solution} = (10 + 20) \text{ cm}^3 = 30 \text{ cm}^3$$

$$\text{Molar concentration of } Cl^- = \frac{1.000 \times 10^{-3}}{30} \times 1000 \text{ M}$$

$$\text{or } [Cl^-] = 3.333 \times 10^{-2} \text{ M}$$

(iv) Calculation of  $K_{sp}$

Now using Eq. 9.3, we can calculate the solubility product of AgCl as shown below:

$$K_{sp} = [Ag^+][Cl^-] = (10^{-6} \times 3.333 \times 10^{-2})$$

$$= \dots\dots\dots$$

9.3.6 Result

The solubility product of silver chloride at ..... °C = .....

Literature value at 25°C =  $1.5 \times 10^{-10}$

9.4 SUMMARY

In this unit, you had studied two experiments – one for determining the pH of solutions and the other for determining the solubility product of a sparingly soluble salt silver chloride. These experiments are based upon the applications of measurement of emf of galvanic cells.



**Hydrochloric acid :  $1.00 \times 10^{-2}$  M :** Using a burette or a measuring cylinder, transfer in small lots  $9 \text{ cm}^3$  of concentrated hydrochloric acid to a  $500 \text{ cm}^3$  volumetric flask which contains about  $250 \text{ cm}^3$  of deionised water. Shake the flask after the addition of each lot. Make up the solution to the mark with deionised water. Let this solution, say solution (A), is approximately 0.1 M. Now titrate this solution against  $1.000 \times 10^{-1}$  M sodium carbonate solution using methyl orange indicator. Calculate the molarity of HCl solution (A) by using the following formula:

$$2 M_{\text{Na}_2\text{CO}_3} \cdot V_{\text{Na}_2\text{CO}_3} = M_{\text{HCl}} \cdot V_{\text{HCl}}$$

Let the molarity of the HCl solution (A) be  $x \text{ mol dm}^{-3}$ .

Prepare HCl solution of  $1.00 \times 10^{-2}$  M concentration by diluting appropriate volume of HCl solution (A) with deionised water. Volume of solution (A) required can be calculated by using the formula,

$$M_1 V_1 = M_2 V_2$$

$$\text{or } V_1 = \frac{M_2 V_2}{M_1}$$

Thus, volume of HCl solution (A) required to prepare  $1 \text{ dm}^3$  of  $1.00 \times 10^{-2}$  M HCl

$$= \frac{1.00 \times 10^{-2} \text{ M} \times 1000 \text{ cm}^3}{x \text{ M}} = \frac{10}{x} \text{ cm}^3$$

Therefore, take  $10/x \text{ cm}^3$  of HCl solution (A) and dilute to  $1 \text{ dm}^3$  with deionised water in a volumetric flask. The concentration of this HCl solution would be  $1.00 \times 10^{-2}$  M.

**Acetic acid :  $1.00 \times 10^{-2}$  M :** Using a burette, transfer  $1.5 \text{ cm}^3$  of glacial acetic acid into a  $250 \text{ cm}^3$  volumetric flask containing  $100 \text{ cm}^3$  of deionised water and make up the solution to the mark. This solution say, acetic acid (A), is approximately 0.1 M. Find out the exact molarity of this solution by titrating against a standard solution of NaOH using phenolphthalein indicator. Let the exact molarity of acetic acid (A) be  $y \text{ mol dm}^{-3}$ .

Now prepare  $1.00 \times 10^{-2}$  M acetic acid by diluting appropriate volume of acetic acid (A) with deionised water. Volume of acetic acid (A) required can be calculated by using molarity equation  $M_1 V_1 = M_2 V_2$ . Thus, volume of acetic acid (A) required to prepare

$$1 \text{ dm}^3 \text{ of } 1.00 \times 10^{-2} \text{ M acetic acid} = \frac{1.00 \times 10^{-2} \text{ M} \times 1000 \text{ cm}^3}{y \text{ M}}$$

$$= \frac{10 \text{ cm}^3}{y}$$

Therefore, take  $10/y \text{ cm}^3$  of acetic acid (A) and dilute with deionised water to  $1 \text{ dm}^3$  in a volumetric flask. This would give an acetic acid solution of concentration  $1.00 \times 10^{-2}$  M.

$1.000 \times 10^{-1}$  M sodium carbonate solution can be prepared by dissolving 1.325 g of anhydrous sodium carbonate in minimum quantity of deionised water and making up to the mark in a  $250 \text{ cm}^3$  volumetric flask.

Details of preparation of standard NaOH solution are given in Appendix-II to Unit 12 of this Block.

**Copper sulphate :  $1.000 \times 10^{-1}$  M** : Dissolve 12.484 g of A.R.  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in deionised water and make up the volume to  $500 \text{ cm}^3$  in a volumetric flask. Before making up the volume, add  $5 \text{ cm}^3$  of acetic acid to prevent hydrolysis.

**Potassium chloride :  $1.000 \times 10^{-1}$  M** : Dissolve 3.728 g of A.R. KCl in deionised water and make up the volume to  $500 \text{ cm}^3$  in a volumetric flask.

**Silver nitrate :  $1.000 \times 10^{-1}$  M** : Dissolve 16.987 g of A.R.  $\text{AgNO}_3$  in deionised water and make up the volume to  $1 \text{ dm}^3$  in a volumetric flask.

**Silver nitrate :  $5.000 \times 10^{-2}$  M** : Dilute  $250 \text{ cm}^3$  of  $1.000 \times 10^{-1}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.

**Silver nitrate :  $1.000 \times 10^{-2}$  M** : Dilute  $50 \text{ cm}^3$  of  $1.000 \times 10^{-1}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.

**Silver nitrate :  $5.000 \times 10^{-3}$  M** : Dilute  $25 \text{ cm}^3$  of  $1.000 \times 10^{-1}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.

**Silver nitrate :  $1.000 \times 10^{-3}$  M** : Dilute  $10 \text{ cm}^3$  of  $5.000 \times 10^{-2}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.

**Silver nitrate :  $5.000 \times 10^{-4}$  M** : Dilute  $25 \text{ cm}^3$  of  $1.000 \times 10^{-2}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.

**Silver nitrate :  $1.000 \times 10^{-4}$  M** : Dilute  $10 \text{ cm}^3$  of  $5.000 \times 10^{-3}$  M  $\text{AgNO}_3$  to  $500 \text{ cm}^3$  with deionised water in a volumetric flask.